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The mastery of the water in soils, its content, its flows in all directions, are so much parameters that it becomes imperative to know with the greatest accuracy. Indeed, when one knows that all polluting rejects are going to finish, sooner or later, by depositing on the ground and penetrate there, helped in that by precipitations, it is essential to know where these waters polluted go. Similarly and so as to master perfectly management of products fertilizing, the conduct of irrigations, etc..., it is fundamental to have the better knowledge as for certain parameter of ground, so;

The hydrodynamics properties of these soils.

*The WIND method allows to determine simultaneously these **hydrodynamics properties** and that on a **SAMPLE** of **GROUND**, so;*

-The relationship "Matric potential-water content $h(\theta)$ "


-The relationship "hydraulic conductivity -water content $k(\theta)$, or conductivity hydraulic-Matric potential $k(h)$."

These fundamental properties govern water transfers in soils.

The WIND method is a LABORATORY method. It consists in submitting a cylinder of soil (collected from the field) initially saturated, to a progressive evaporation and to follow:

-The changes in the time of the mass of soil $m(t)$ in the cylinder so as to determine the water loss by evaporation and the average humidity of the ground.

-the evolution of the value of matric potential $h(z,t)$ soil, at different depths (z) in function of the time.

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
CENTRALIZING BOX

This device (centralizing box) contain the differential pressure sensors (x6) and the electronic circuit. Of this "electronic box" come out six micro-tensiometers to put it in the Stainless Steel sample ring which contain the soil to analyze. This sample RING (which contain soil) is put on an very accurate balance. This balance will allow to know the water content inside the sample ring with a great precision at any time. The knowledge of water content and hydric potentials (given by micro-tensiometers) in the same time inside the ring, allow to establish the relationship between both. The centralizing box was designed for to lessen all disturbances as external temperature, external pressure etc...

This device allows to do **very accurate measurements** on undisturbed soil **from the field** but in the **laboratory**.

SYSTEME "WIND" (Caractérisation hydrodynamique des sols)**"WIND" system (Hydrodynamic characterization of soils)****WIND****SPECIFICATIONS ELECTRIQUES**

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ELECTRICAL SPECIFICATIONS

- **Power supply** : 12 volts c.c (9 to 18 V possible.).
- **Output voltage for pressure sensor** : 0 to 50 mV.
(These values can be different for each sensor. The calibration certificate given with each points the measurement range for each sensors with a great accuracy, and the correspondence hPa (mbar).
- **Voltage regulator for pressure sensors** : 5,00 Volts / 400 mA.
Voltage output : $\pm 2\%$.
Line regulation : $\pm 0,5\%$ maxi.
Charge regulation : $\pm 1\%$ maxi. (10-90%).
Residual ripples : (20Mhz pass-band) : 60 mV.
typical crest to crest.
Protection against short-circuits : continuous.
Operating temperature : -10°C to +40°C. (ATTENTION, the centralizing box does not to be at these temperatures).
- **Temperature of centralizing box** : from +15°C to + 25°C.
- **Internal measurement of box temperature** : by PT100 class A.
- **Inner resistances network** : Welwyn type 15 ppm / 0,1%.
- **Differential pressure sensor** :
Offset : 0 $\pm 1,5$ mV (given by calibration certificate).
Linearity : $\pm 0,25$ (%) SPAN).
Repeatability & hysteresis : $\pm 0,20$ (%) SPAN).
Response time : 1 msec.
Stability over one year : $\pm 0,5$ (%) SPAN).
Compensated temperature : from 0°C to +50°C.
Sensitivity : 6,67mV/psi (0,1mV/mbar).
Overpressure input : + 2,5 bars.
- **Precision balance** : Consult us.

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DESCRIPTION of SYSTEM

- **Centralizing box** : (reference : 7703000)
Price : see page 88
This device consists of :
 - One watertight box with transparent lid (225 x 175 x 90 mm).
 - One electronic card for the sensors framing.
 - Six differential pressure sensors (-1; 0 ;+1 bar). Compensation of external temperature.
 - Six ceramic pipes (SDEC220 type (*)) sticked on special capillary (\varnothing : 3 x 1,5 mm) with all silicon coupling sleeves.
 - One external RS 232 type plug for connecting output of the centralizing box to the datalogger (datalogger non include).
 - One calibration certificate for each sensor.
- **Precision balance** : consult us.
- **Software (Window 95-98)** : In progress

Field EQUIPMENTS

- **Soil sampling cylinder: set of 3 pieces** (réf : 7710114).
(ATTENTION : this price include 3 cylinders and their lids).
- **Cylinder holder with hammering head**: ref 7700115
- **Spirit level** : Consult us.
- **Club hammer impact absorbing type**: ref 040505